

# EasyTraceCleaner® (ETC)

*Designed to efficiently clean your containers, bottles, microwave liners and pipette tips for trace or ultratrace analysis*



- ✓ Acid resistant
- ✓ Optimal efficiency
- ✓ Economical
- ✓ Easy to use
- ✓ Safe, use 24/7
- ✓ Maintenance free



Specialist Chemicals & Products for Chemistry

[sales@qmx.com](mailto:sales@qmx.com) [www.qmx.com](http://www.qmx.com)


Qmx Laboratories, Thaxted, CM6 2PY T 01371 831611

# EasyTraceCleaner<sup>®</sup>

- From 30ml to 250ml PFA-PP vials and jars
- Bottles
- Microwave liners

Up to 48 vials + caps per cleaning cycle

- ✓ Low acid consumption
- ✓ Standard grade acid only used
- ✓ Offered in different formats to suit all vessel sizes



Optional plate to double your cleaning batches (up to 100)



*Temperature controller must be ordered separately*



Specialist Chemicals & Products for Chemistry

sales@qmx.com www.qmx.com

Qmx Laboratories, Thaxted, CM6 2PY T 01371 831611