

# PremiSolv™ - ICP Solvent

***“The Smell of Success - Is No Smell”  
- No Sulfur/Aromatic Odour***



PremiSolv™ ICP solvent is a replacement for low odour kerosene and xylene for use as a diluent in ICP/DCP analysis of metals in oil and other organic fluids. PremiSolv™ ICP solvent is supplied in one gallon containers for shipment via common carrier or five gallons for high volume users. PremiSolv™ ICP solvent features:

- ❖ Supplied with certificate of analysis listing the concentration of 34 different metals and sulphur
- ❖ Extremely low toxicity compared to kerosene or xylene
- ❖ No sulfur/aromatic odour for increased worker comfort and less exposure
- ❖ Non-Hazardous for shipping/disposal\*
- ❖ Dilution of samples and standards

\* flash point >200°F / 93.33°C



***Free sample  
visit [conostan.com](http://conostan.com)  
for more details***

Description	Catalogue Number
PremiSolv™, 5 Gallons	150-700-002
PremiSolv™, 1 Gallon	150-700-003
PremiSolv™, Sample 400 mL	150-700-000

**CONOSTAN®**  
Oil Analysis Standards

***“The benchmark for  
accurate measurements”***

qm<sub>x</sub>

+44 (0)1371 831611  
info@qmx.com  
www.qmx.com